

BODYVIZ

Hardware and Software Requirements

July 1, 2009

*** CONFIDENTIAL ***

Hardware Requirements and Recommendation

Item	Recommended Specification	Minimum Specification
Processor	Intel Core 2 Duo or similar	Pentium 4, 1.0 GHz or above
Hard disk	300MB of free hard disk space for software and supporting libraries. Additional hard disk space will be needed for imaging data.	300MB of free hard disk space for software and supporting libraries. Additional hard disk space will be needed for imaging data.
RAM	2 GB or better if available	1 GB
Graphics	A video graphics card (i.e. NVIDIA) with 512MB or more of video RAM	A video graphics card (i.e. NVIDIA) with 256MB of video RAM (128MB will run but performance will suffer)
Controller Unit	Microsoft Xbox 360 Wired or Wireless Controller with USB PC adapter (Use of wireless depends on your facility's RF emissions)	Microsoft Xbox 360 Wired or Wireless Controller with USB PC adapter (Use of wireless depends on your facility's RF emissions)

Note: Computer hardware safety and emissions requirements are specified by your facility.

Software Requirements and Recommendation

Item	Recommended Specification	Minimum Specification
Operating System	Windows XP, or Vista	Windows 2000 with SP4
Software	Microsoft .NET v3.5 (If Microsoft .NET 3.5 is not already on computer, BodyViz installer will download and install necessary components from Microsoft.)	Microsoft .NET v3.5 (If Microsoft .NET 3.5 is not already on computer, BodyViz installer will download and install necessary components from Microsoft.)